



The future of vision

SIOBAL



SIOBAL Silicone Oil

Reference	S1000	S2000	S5000
Brand name	SIOBAL SILICONE OIL		
Definition	Buffer agent for retinal detachment.		
Material	High-purity silicone oil free of low-molecular-weight components.		
Viscosity	1 000 cSt	2 000 cSt	5 000 cSt
Indications	<ul style="list-style-type: none"> • Retinal detachment. • Proliferative vitreoretinopathy. • Perforating wounds. • Proliferative diabetic retinopathy. 		
Sterilization method	Steam		
Supply	Sterile, single unit. 10ml prefilled glass syringe.		
Shelf-life	3 years		

CE 0318

REV: 01